

# TECNIS<sup>™</sup> Toric Calculator

## User Tutorial

Welcome to the TECNIS<sup>™</sup> Toric Calculator.

The TECNIS<sup>™</sup> Toric Calculator has been specifically designed to make calculating the TECNIS<sup>™</sup> Acrylic 1-Piece Toric Intraocular Lenses (IOLs) accurate, intuitive, and easy.

With the TECNIS<sup>™</sup> Toric Calculator you can:

- Input prompted data in order to determine the TECNIS<sup>™</sup> Toric IOL model that will provide an optimized surgical outcome.
- Get calculated TECNIS<sup>™</sup> Toric IOL results for your surgical needs.
- Print a report of patient information and results for surgeons' reference during cataract surgery.

# TECNIS™ Toric Calculator

## Overview

Johnson & Johnson Surgical Vision, Inc. (Johnson & Johnson Vision) has developed the Toric Intraocular Lens (IOL) models with patients undergoing cataract surgery in mind. The goal is to deliver a lens which will reduce residual refractive astigmatism based on a personalized level of data entry. The TECNIS™ Toric Calculator considers preoperative corneal astigmatism, as well as any anticipated surgically induced astigmatism in determining the optimal Toric IOL power for the correction of postoperative corneal astigmatism. The TECNIS™ Toric Calculator also provides the option of including an algorithm to account for posterior corneal astigmatism. The practice that will afford appropriate results requires that accurate keratometry measurements are used to determine preoperative corneal astigmatism. While the desired surgical procedure strives towards astigmatic neutrality, even minor incisions during cataract surgery can induce corneal astigmatism. This may affect the amount of postoperative corneal astigmatism as well as the respective axis. As a result, Johnson & Johnson Vision recommends that based upon individual surgical technique and past results, surgeons should customize their surgically induced corneal astigmatism value.

The TECNIS™ Toric Calculator is a tool, designed to help further integrate toric lenses into your practice. It is a software application that helps you select an appropriate TECNIS™ IOL for your patient. In the TECNIS™ Toric Calculator, the spherical equivalent IOL power is introduced as an input for each patient. Surgeons are able to use their spherical equivalent IOL power calculation method of choice. The TECNIS™ Toric Calculator calculates different cylinder IOL power options for you, as well as the orientation in which the IOL should be implanted. In addition, predicted postoperative residual astigmatism is calculated for each cylinder IOL power suggested and, placed in the indicated orientation in the patient's eye in order to facilitate the surgeon's lens selection.

In the TECNIS™ Toric Calculator, the toricity and residual astigmatism are calculated for each particular eye in a customized fashion, using the Holladay 1 formula instead of using a fixed ratio based on average ocular physiology. The TECNIS™ Toric Calculator allows for the input of personalized surgically induced astigmatism and/or IOL A-constants, as well as different options for keratometry input, to fit the needs of your practice.

Although the tool has been tested to a large range of possible situations, it is the user's responsibility to decide the usefulness of the calculated information.

## Intended Use / Purpose

The TECNIS™ Toric Calculator is a software application that informs the selection of an appropriate TECNIS™ IOL.

## Indications for Use

The TECNIS™ Toric Calculator is a software application that informs the selection of an appropriate TECNIS™ Toric IOL that is indicated to correct aphakia following removal of the natural crystalline lens and treatment of preexisting astigmatism.

## Intended Users

The intended users for the TECNIS™ Toric Calculator are clinicians, ophthalmologists, and/or eye surgeons.

## Contraindications

None.

## Expected Clinical Benefits

The TECNIS™ Toric Calculator assists surgeons in selecting the appropriate TECNIS™ IOL Model to fit patient needs.

# TECNIS™ Toric Calculator

## Prerequisites

To use the TECNIS™ Toric Calculator, you will need: An Apple® computer (Mac OS® 14) or PC (Windows® 10 operating system).

- A working Internet connection
- Safari®, Google Chrome™ or Microsoft Edge® browser

## Documentation Prerequisites

- Adobe® Reader® Software ([Download it here.](#))

**Note:** This manual is applicable for all approved TECNIS™ Toric IOLs and their corresponding non-Toric IOLs; TECNIS™ IOL is listed for example throughout this manual.

## Performance Characteristics

The performance and safety of the TECNIS™ Toric Calculator has been established through pre-clinical testing and clinical use and represents state-of-the-art medical software. The TECNIS™ Toric Calculator performs as intended for the function as a software application that helps to calculate and select the appropriate TECNIS™ IOL for postoperative corneal astigmatism correction based on user inputs when used according to the instructions for use and labeling.

## Warnings & Precautions

None.

Note: Warnings and precautions for the specific intraocular lenses are described in their respective instructions for use.

Screen images shown below are for illustration purposes only.

# Quick Setup Guide

Read this page if this is your first time using the TECNIS™ Toric Calculator.

Within this guide, you will find the recommended first steps for getting your TECNIS™ Toric Calculator configured and ready to perform calculations.

After accessing the TECNIS™ Toric Calculator for the first time, a Registration page, similar to the one shown highlighted to the right, is displayed.

Once you are on the Registration page:

1. Select your location and language from the drop-down lists.
2. Review and Accept the TECNIS™ Toric Calculator Terms of Service.

Once the Terms of Service are accepted the Calculator page is displayed, and the TECNIS™ Toric Calculator is ready for use.

## Registration Page

TECNIS™ Toric Calculator

TECNIS™ Family of IOLs

Johnson & Johnson VISION

Registration

Europe - Middle East - Africa English I have read and agree to the terms below.

Accept

TECNIS™ intraocular lenses are currently available for sale in select markets. This application is therefore intended only for residents of those markets in which the respective TECNIS™ intraocular lens(es) are available, and this application may only be used for those TECNIS™ intraocular lens(es) that are approved for use in your market. Using this mobile application and clicking the "Accept" button constitutes your acknowledgement of these restrictions, and your acknowledgement that you selected your proper location of residence.

TECNIS™ TORIC CALCULATOR LICENSE AGREEMENT, TERMS OF USE AND PRIVACY POLICY

PLEASE READ THIS USE AND LICENSE AGREEMENT ("AGREEMENT") IN FULL

TERMS OF SERVICE LRI CALCULATOR

The product-specific application is intended for the residents of a particular country or countries. As a result, the application may contain information on pharmaceuticals, medical devices and other products or uses of those products that are not approved in other countries or regions. If you are a resident of a country other than those to which the application lists, please return to [surgical.jnjvision.com](http://surgical.jnjvision.com) or contact your local Johnson & Johnson Vision affiliate to obtain the appropriate product information for your country of residence.

This application is intended for use in connection with Johnson & Johnson Vision IOLs only and should not be used in connection with any other IOLs.

©2024 Johnson & Johnson Surgical Vision, Inc.

Unless otherwise specified, all product names appearing in this program are trademarks owned by or licensed to Johnson & Johnson Vision, its subsidiaries or affiliates. No use of any Johnson & Johnson Vision trademark, trade name, or trade dress in this program may be made without the prior written authorization of Johnson & Johnson Vision, except to identify the product or service of the company.

# Lens Selection

If there are more than one type of TECNIS™ Toric IOL available for the selected location, the TECNIS™ Family of IOLs logo is displayed in the upper-left corner until a lens is selected. If there is only a single TECNIS™ Toric IOL available for the selected location, the lens will be automatically selected and the logo for the lens is displayed in the upper-left corner.

To select from more than one TECNIS™ Toric IOL, the “SELECT LENS” button displays the available lenses for the current location. Simply select the lens of choice from the list of lenses displayed.

The logo for the newly selected lens will be displayed in the upper-left corner, and the “SELECT LENS” button is now labeled “CHANGE LENS” but functions the same.

## Calculator Page – Lens Selection

The screenshot displays the TECNIS™ Toric Calculator interface. The top left corner features the TECNIS™ Family of IOLs logo. The main area is divided into several sections:

- Surgeon and Patient Information:** Includes fields for Surgeon Name, Patient Information, Eye Selection (OD/Right or OS/Left), and K Notation.
- Keratometry:** Includes fields for Surgically Induced Astigmatism (SIA), Flat K1, Steep K2, and Preop Corneal Astigmatism.
- Biometry:** Includes fields for Axial Length, Method (Optical or Immersion), and A-constant.
- Calculator:** A central panel showing a list of available lenses:
  - TECNIS Eyhance™ with TECNIS SIMPLICITY™ (DIU)
  - TECNIS™ Toric II 1-Piece IOL (ZCU)
  - TECNIS™ Toric 1-Piece IOL (ZCT)
  - TECNIS PureSee™ with TECNIS SIMPLICITY™ (DET)
  - TECNIS Synergy™ with TECNIS SIMPLICITY™ (DFW)
  - TECNIS Symphony™ with TECNIS SIMPLICITY™ (DXW)
  - TECNIS Symphony™ (ZXT)
  - TECNIS Multifocal (ZMT)
- Results:** A table showing IOL Details and Residual Astigmatism.

The interface also includes a circular diagram for astigmatism correction with angles (0°, 45°, 90°, 135°, 180°, 225°, 270°, 315°) and a "N a s a l" label. The bottom of the page contains links for "TERMS OF SERVICE" and "LRI CALCULATOR".

This is a representative image of lens selection and not all lens options are available worldwide.

# Calculation Results

Data entry fields, shown highlighted to the right, are on the left side of the page. However, the Eye Selection can also be made by selecting the appropriate eye on the Patient Face on the right side of the page.

Select any of the Section Headers with the circled information symbols to display definition and data input range information for each field within that section.

**Note:** Make sure to always input anterior corneal astigmatism measurements, regardless of “Include Posterior Corneal Astigmatism” checkbox selection. To ensure accurate calculation results, total corneal astigmatism values should never be used.

**Note:** A predetermined value for posterior corneal astigmatism can be included in the calculation by checking the box labeled “Include Posterior Corneal Astigmatism”. The option to include the predetermined value of posterior corneal astigmatism is based on an algorithm that combines published literature<sup>1</sup> with a retrospective analysis of existing clinical data.

**Note:** Personalized constants for the corresponding non-Toric model (i.e., Model ZCB00) are valid starting points for Toric model calculations. Note that the formula used for cylinder power calculation is the Holladay-1.

Once the data entry is complete, select the “Calculate Results” button to calculate the results.

## Calculator Page – Data Entry

The screenshot displays the TECNIS™ Toric Calculator interface. At the top, it shows the TECNIS Eyhance™ logo, the text 'DIU', and 'Europe - Middle East - Africa'. A 'CHANGE LENS' dropdown menu is visible. The interface is divided into two main sections: 'Calculator' (left) and 'Preferences & Support' (right).

The 'Calculator' section contains several data entry fields:

- Surgeon and Patient Information:** Surgeon Name, Date (11/14/2024), Patient Information, Patient Age, Eye Selection (OD (Right) selected), OS (Left), K Notation (D selected), and mm.
- Keratometry:** Surgically Induced Astigmatism (SIA) in D, @ Meridian (Incision Location), Flat K1 in mm, Flat K1 @ Meridian, Steep K2 in mm, Steep K2 @ Meridian, Preop Corneal Astigmatism in D, and an 'Include Posterior Corneal Astigmatism' checkbox.
- Biometry:** Axial Length in mm, Method (Optical or Immersion), and A-constant (119.30).
- Calculation Preferences:** SE IOL Power (-Select-), K Index (-Select-), and Refractive Cylinder Convention (Plus selected, Minus).

The 'Preferences & Support' section features a patient face diagram for eye selection (OD highlighted), a circular diagram for astigmatism orientation with meridians at 90°, 45°, 0°, 315°, 270°, 225°, 180°, and 135°, and a 'Final Results' table.

The 'Final Results' table is structured as follows:

IOL Details		Residual Astigmatism	
IOL Model	Orientation	Cylinder	Axis

At the bottom of the calculator section, there are buttons for 'Calculate Results', 'Clear Entries', and 'Print Results'. The footer contains 'TERMS OF SERVICE' and 'LRI CALCULATOR'.

<sup>1</sup> Koch DD et al. Contribution of posterior corneal astigmatism to total corneal astigmatism. J Cataract Refract Surg. 2012 Dec;38(12):2080-7

# Calculation Results

Calculation results, shown highlighted to the right, are displayed on the right side of the page.

The postoperative corneal astigmatism orientation (orange / red: steep axis) and incision location (blue) are displayed on the diagram of the eye located in the graphical area on the right side of the page, as well as the recommended IOL orientation. Identification of Temporal and Nasal are displayed, relative to the choice of eye (OD or OS).

The Final Results table contains different Toric IOL cylinder power suggestions for this eye with their respective orientations and predicted residual astigmatism, in the convention chosen by the Calculation Preferences section.

Select the Final Results section header for more information.

Once an IOL Model is selected from the Final Results, the "Print Results" button is enabled to allow printing of the results.

## Calculator Page – Calculation Results

**TECNIS Eyhance™** DIU  
with TECNIS SIMPLICITY™

TECNIS™ Toric Calculator  
Europe - Middle East - Africa  
CHANGE LENS ▾

Johnson & Johnson VISION

Calculator Preferences & Support

**Surgeon and Patient Information** Notes...

Surgeon Name: Surgeon Date: 11/14/2024

Patient Information: XYZ-123-ABC Patient Age: 65

Eye Selection:  OD (Right)  OS (Left) K Notation:  D  mm

**Keratometry**

Surgically Induced Astigmatism (SIA): 0.30 D @ Meridian (Incision Location): 180

Flat K1: 41.25 D Flat K1 @ Meridian: 12

Steep K2: 43.34 D Steep K2 @ Meridian: 102

Preop Corneal Astigmatism: 2.09 D  Include Posterior Corneal Astigmatism

**Biometry**

Axial Length: 24.11 mm

Method: Optical or Immersion

A-constant: 119.30

**Calculation Preferences**

SE IOL Power: 20.0 D

K Index: 1.3375

Refractive Cylinder Convention:  Plus  Minus

Calculate Results Clear Entries Print Results

**Final Results**

IOL Details		Residual Astigmatism	
IOL Model	Orientation	Cylinder	Axis
<input type="radio"/> DIU300	101°	+0.33 D	101°
<input type="radio"/> DIU375	101°	+0.16 D	11°
<input type="radio"/> DIU450	101°	+0.65 D	11°

TERMS OF SERVICE LRI CALCULATOR

# Print Results

**Note:** User-entered data and calculated results are not stored when the user exits the TECNIS™ Toric Calculator.

In order to keep a record of the TECNIS™ Toric Calculator output for the patient's file and for surgeon reference during the procedure:

1. Select an IOL Model from the Final Results table on the Calculator page.
2. Select the "Print Results" button at the bottom of the Calculator page.
3. Verify the input and results and proceed according to the confirmation message displayed.
4. When the Print dialog is displayed, select the printer of choice, and select the "OK" button to print and display the Results page.
5. Select the "Return to Calculator" button to leave the Results page and return to the Calculator page.

## Print Results Page




DIU

[Return to Calculator](#)



Surgeon Name	Patient Information	Date	Patient Age
Surgeon	XYZ-123-ABC	2/6/2024	65

OD (Right)





— (Incision Location) 180°

— Calculated Orientation

IOL Details		Residual Refraction	
IOL Model	Orientation	Cylinder	Axis
DIU300	101 °	+0.33 D	101 °

Calculation Details			
SE IOL Power	K Index	Refractive Cylinder Convention	
+20.0	1.3375	Plus	

Keratometry and Biometry						
Flat K1 @ Meridian	Steep K2	SIA @ Meridian	AxL Length	Biometry Method	A-constant	PCA
41.25@12°	43.34@102°	0.30@180°	24.11	Optical or Immersion	119.30	No

IOL Details		Residual Astigmatism	
IOL Model	Orientation	Cylinder	Axis
DIU300	101 °	+0.33 D	101 °
DIU375	101 °	+0.16 D	11 °
DIU450	101 °	+0.65 D	11 °

Notes

MD CE UDI (01)05050474708822(8012) - 4.9.0.0

# Clearing Data Input and Calculation Results

If you would like to initiate a new calculation, select the “Clear Entries” button to clear all previously entered data and calculation results. A prompt to confirm your choice, shown highlighted to the right, is displayed.

**Note:** Changing the Patient Information or Patient Age without clearing previous data and results will display a similar prompt to confirm whether you want to clear all the data entries or not.

**Note:** Changing the selected Lens Type without clearing previous data and results will also display a similar prompt to confirm whether you want to clear all the data entries or not.

It is good practice to choose to clear all data and results when entering information for a new patient.

## Calculator Page – Clear Entries

The screenshot shows the LRI Calculator interface with a 'Clear Entries' dialog box open. The dialog box contains the following text: 'Clear Entries', 'Please confirm that you wish to clear all entries that you have made on this screen. (e.g., OD (Right) & OS (Left))', and 'OK' and 'Cancel' buttons. The background interface includes the following sections:

- Surgeon and Patient Information:** Surgeon Name (Surgeon), Date (11/14/2024), Patient Information (XYZ-123-ABC), Patient Age (65), Eye Selection (OD (Right) selected, OS (Left) unselected), K Notation (D selected, mm unselected).
- Keratometry:** Surgically Induced Astigmatism (SIA) (0.30 D @ Meridian (Incision Location) 180), Flat K1 (41.25 D Flat K1 @ Meridian 12), Steep K2 (43.34 D Steep K2 @ Meridian 102), Preop Corneal Astigmatism (2.09 D), Include Posterior Corneal Astigmatism (checkbox unselected).
- Biometry:** Axial Length (24.11 mm), Method (Optical or Immersion), A-constant (119.30).
- Calculation Preferences:** SE IOL Power (20.0 D), K Index (1.3375), Refractive Cylinder Convention (Plus selected, Minus unselected).
- Final Results:** IOL Details table with columns: IOL Model, Orientation, Residual Astigmatism (Cylinder, Axis).

At the bottom of the interface, there are buttons for 'Calculate Results', 'Clear Entries', and 'Print Results'. The footer contains 'TERMS OF SERVICE' and 'LRI CALCULATOR'.

# Preferences and Support - Preferences

Preferences are accessed by selecting the “Preferences & Support” button. The Preferences & Support page has two sections. The upper section, shown highlighted to the right, is for setting the user preferences.

The Preferences are separated into settings for the Display and for the Calculation.

The Display preferences allow the user to change the default location and language. In addition, the K Notation can be set to default to either “mm” for millimeters or “D” for diopters and the Decimal format can be set to default to “0.00” or “0,00” for numeric entries.

The Calculation preferences allow the user to change the default settings for the Refractive Cylinder Convention, the K Index and whether to Include Posterior Corneal Astigmatism. In addition, the default Biometry Method can be selected. If a personalized constant is selected for the Biometry Method the corresponding value field allows the user to enter their personalized A-Constant, Surgeon Factor or Anterior Chamber Depth accordingly.

Select the “Apply Preferences” button, on the right side of the Preferences header, to save any changes made to the preferences.

Select the “Reset to Defaults” button, on the right side of the Preferences header, to reset all of the preferences, except Location and Language, to the system defaults.

## Preferences and Support – Preferences Page

**TECNIS™** Toric Calculator  
Europe - Middle East - Africa

Johnson & Johnson VISION

Calculator | Preferences & Support

**Preferences** Reset To Defaults Apply Preferences

**Display**

Location: Europe - Middle East - Africa

Language: English

K Notation:  D  mm

Decimal:  0.00  0,00

**Calculation**

Refractive Cylinder Convention:  Plus  Minus

K Index: -Select-

Method: Optical or Immersion

A-constant: 119.30

Include Posterior Corneal Astigmatism

**Support**

**Training and Information Documents**

- [TECNIS™ Toric Calculator User Tutorial](#)
- [TECNIS™ Toric Calculator Frequently Asked Questions \(FAQ\)](#)
- [TECNIS™ Toric Calculator Terms of Service](#)

www.e-ifu.com EU: +353 1 7782500

Johnson & Johnson Surgical Vision, Inc.  
31 Technology Drive, Suite 200  
Irvine, CA 92618 USA  
2024-11-14

AMO Ireland  
Block B  
Liffey Valley Office Campus  
Quarryvale, Co. Dublin, Ireland

MD CE 0344 UDI (01)05050474708822(8012) - 4.9.0.0

TERMS OF SERVICE | LRI CALCULATOR

# Preferences and Support - Support

Support documents are accessed by selecting the "Preferences & Support" button. The Preferences & Support page has two sections. The lower section, shown highlighted to the right, is for accessing the support documents, and indicating compliance information along with the version number of the application.

There are three linked documents in the Support section. Please refer to the Documentation Prerequisites shown previously for a link to install the required Adobe® Reader® software if it is not already installed.

Note: The Version Number uniquely identifies the software on the site and is displayed for reference in the lower-right corner of the Support section. The date shown below the manufacturer address in the upper-right corner of the Support section corresponds to the date the software version was built.

## Preferences and Support – Support Page

**TECNIS™** Family of IOLs  
TECNIS™ Toric Calculator  
Europe - Middle East - Africa  
Johnson & Johnson VISION

Calculator | Preferences & Support

Preferences Reset To Defaults Apply Preferences

**Display**  
Location: Europe - Middle East - Africa | Language: English | K Notation:  D  mm | Decimal:  0.00  0,00

**Calculation**  
Refractive Cylinder Convention:  Plus  Minus | K Index: -Select- | Method: Optical or Immersion | A-constant: 119.30  
 Include Posterior Corneal Astigmatism

**Support**

**Training and Information Documents**  
[TECNIS™ Toric Calculator User Tutorial](#)  
[TECNIS™ Toric Calculator Frequently Asked Questions \(FAQ\)](#)  
[TECNIS™ Toric Calculator Terms of Service](#)

Johnson & Johnson Surgical Vision, Inc.  
31 Technology Drive, Suite 200  
Irvine, CA 92618 USA  
2024-11-14

AMO Ireland  
Block B  
Liffey Valley Office Campus  
Quarryvale, Co. Dublin, Ireland

www.e-ifu.com EU: +353 1 7782500

TERMS OF SERVICE | LRI CALCULATOR

## Serious Incidents Reporting

For a patient/user/third party in the European Union and in countries with identical regulatory regime; if, during the use of this software or as a result of its use, a serious incident has occurred, please report it to Johnson & Johnson Surgical Vision, Inc. and/or its authorized representative and to your national authority.

### Symbol/Title (explanation):

SYMBOL	TITLE (EXPLANATION)
	Manufacturer. The date below the manufacturer name and address is the Date of manufacture (YYYY-MM-DD: year-month-day).
	Authorized representative in the European Community/European Union
	Authorized representative in Switzerland
	Consult instructions for use or consult electronic instructions for use
	Medical Device
	Unique Device Identifier
	Importer
	European Community/European Union
	Switzerland



Johnson & Johnson Surgical Vision, Inc.  
31 Technology Drive, Suite 200,  
Irvine, CA 92618 USA

 AMO Ireland  
Block B  
  Liffey Valley Office Campus  
Quarryvale, Co. Dublin, Ireland

 AMO Switzerland GmbH  
Gubelstrasse 34  
  6300 Zug  
Switzerland

  
0344

Supplied in Australia and New Zealand by:

AMO Australia Pty Ltd.  
1-5 Khartoum Road,  
North Ryde, NSW 2113  
Australia

TECNIS is a trademark owned by or licensed to Johnson & Johnson Surgical Vision, Inc.

All other trademarks are the intellectual property of their respective owners.

© 2024 Johnson & Johnson Surgical Vision, Inc.