Calculadora tórica TECNIS[™]

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Johnson & Johnson VISION

Registration

Europe - Middle East - Africa English I have read and agree to the terms below.

Accept

TECNIS™ intraocular lenses are currently available for sale in select markets. This application is therefore intended only for residents of those markets in which the respective TECNIS™ intraocular lens(es) are available, and this application may only be used for those TECNIS™ intraocular lens(es) that are approved for use in your market. Using this mobile application and clicking the "Accept" button constitutes your acknowledgement of these restrictions, and your acknowledgement that you selected your proper location of residence.

TECNIS™ TORIC CALCULATOR LICENSE AGREEMENT, TERMS OF USE AND PRIVACY POLICY

PLEASE READ THIS USE AND LICENSE AGREEMENT ("AGREEMENT") IN FULL

TERMS OF SERVICE LRI CALCULATOR

The product-specific application is intended for the residents of a particular country or countries. As a result, the application may contain information on pharmaceuticals, medical devices and other products or uses of those products that are not approved in other countries or regions. If you are a resident of a country other than those to which the application lists, please return to surgical.jnjvision.com or contact your local Johnson & Johnson Vision affiliate to obtain the appropriate product information for your country of residence.

This application is intended for use in connection with Johnson & Johnson Vision IOLs only and should not be used in connection with any other IOLs.

©2024 Johnson & Johnson Surgical Vision, Inc.

Unless otherwise specified, all product names appearing in this program are trademarks owned by or licensed to Johnson & Johnson Vision, its subsidiaries or affiliates. No use of any Johnson & Johnson Vision trademark, trade name, or trade dress in this program may be made without the prior written authorization of Johnson & Johnson Vision, except to identify the product or service of the company.

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The screenshot displays the TECNIS™ Toric Calculator interface. The top left corner features the TECNIS™ Family of IOLs logo. The main area is divided into several sections: Surgeon and Patient Information, Keratometry, and Biometry. A central panel titled "SELECT LENS" lists various lens options, including TECNIS Eyhance™ DIU, TECNIS Toric II 1-Piece IOL ZCU, TECNIS Toric 1-Piece IOL ZCT, TECNIS PureSee™ DET, TECNIS Synergy™ DFW, TECNIS Symphony™ DXW, TECNIS Symphony™ ZXT, and TECNIS Multifocal ZMT. A circular diagram on the right side shows the orientation of the lens, with angles ranging from 0° to 315°. The bottom right corner displays the Results section, which includes IOL Details and Residual Astigmatism. The interface also includes a "Calculate Results" button and a "Print Results" button.

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The screenshot displays the 'TECNIS™ Toric Calculator' interface. The top navigation bar includes the TECNIS Eyhance logo, 'DIU', 'Europe - Middle East - Africa', and the Johnson & Johnson VISION logo. Below this is a 'CHANGE LENS' dropdown menu. The main interface is divided into two tabs: 'Calculator' (active) and 'Preferences & Support'. The 'Calculator' tab contains several sections: 'Surgeon and Patient Information' with fields for Surgeon Name, Date (11/14/2024), Patient Information, Patient Age, Eye Selection (OD/OS), K Notation (D/mm), and a 'Notes...' button; 'Keratometry' with fields for Surgically Induced Astigmatism (SIA), Flat K1, Steep K2, Preop Corneal Astigmatism, and an 'Include Posterior Corneal Astigmatism' checkbox; 'Biometry' with fields for Axial Length, Method (Optical or Immersion), and A-constant (119.30); and 'Calculation Preferences' with dropdowns for SE IOL Power, K Index, and a 'Refractive Cylinder Convention' (Plus/Minus) radio button. A circular diagram on the right shows meridians from 0° to 315° with 'T e m p o r a l' and 'N a s a l' labels. The 'Final Results' section contains a table for IOL Details and Residual Astigmatism. At the bottom are 'Calculate Results', 'Clear Entries', and 'Print Results' buttons, along with 'TERMS OF SERVICE' and 'LRI CALCULATOR' links.

¹ Koch DD et al. Contribution of posterior corneal astigmatism to total corneal astigmatism. J Cataract Refract Surg. 2012 Dec;38(12):2080-7

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The screenshot displays the TECNIS™ Toric Calculator interface. The top navigation bar includes the TECNIS Eyhance™ logo, the text 'DIU', 'Europe - Middle East - Africa', and the Johnson & Johnson VISION logo. Below this is a 'CHANGE LENS' button. The main interface is divided into two sections: 'Calculator' and 'Preferences & Support'.

Calculator Section:

- Surgeon and Patient Information:** Surgeon Name (Surgeon), Date (11/14/2024), Patient Information (XYZ-123-ABC), Patient Age (65). Eye Selection: OD (Right) selected, OS (Left) unselected. K Notation: D selected, mm unselected.
- Keratometry:** Surgically Induced Astigmatism (SIA) (0.30 D), Flat K1 (41.25 D), Steep K2 (43.34 D), Preop Corneal Astigmatism (2.09 D). Meridian (Incision Location) (180), Flat K1 @ Meridian (12), Steep K2 @ Meridian (102). Include Posterior Corneal Astigmatism (checkbox unselected).
- Biometry:** Axial Length (24.11 mm), Method (Optical or Immersion), A-constant (119.30).
- Calculation Preferences:** SE IOL Power (20.0 D), K Index (1.3375), Refractive Cylinder Convention (Plus selected, Minus unselected).

Preferences & Support Section:

- Diagram of the Eye:** Shows a circular diagram with meridians at 90°, 135°, 180°, 225°, 270°, 315°, and 45°. A blue arc indicates the incision location at 180°. A red arc indicates the astigmatism axis at 101°. Labels 'T e m p o r a l' and 'N a s a l' are present.
- Final Results Table:**

IOL Details		Residual Astigmatism	
IOL Model	Orientation	Cylinder	Axis
<input type="radio"/> DIU300	101 °	+0.33 D	101 °
<input type="radio"/> DIU375	101 °	+0.16 D	11 °
<input type="radio"/> DIU450	101 °	+0.65 D	11 °

Buttons at the bottom: Calculate Results, Clear Entries, Print Results. Footer: TERMS OF SERVICE, LRI CALCULATOR.

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DIU


[Return to Calculator](#)

Surgeon Name	Patient Information	Date	Patient Age
Surgeon	XYZ-123-ABC	2/6/2024	65

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The screenshot displays the LRI Calculator interface with a 'Clear Entries' dialog box overlaid. The dialog box contains the following text: 'Clear Entries', 'Please confirm that you wish to clear all entries that you have made on this screen. (e.g., OD (Right) & OS (Left))', and 'OK' and 'Cancel' buttons. The background interface includes the following sections:

- Surgeon and Patient Information:** Surgeon Name (Surgeon), Date (11/14/2024), Patient Information (XYZ-123-ABC), Patient Age (65).
- Eye Selection:** OD (Right) selected, OS (Left) unselected. K Notation: D selected, mm unselected.
- Keratometry:** Surgically Induced Astigmatism (SIA) 0.30 D @ Meridian (Incision Location) 180. Flat K1 41.25 D Flat K1 @ Meridian 12. Steep K2 43.34 D Steep K2 @ Meridian 102. Preop Corneal Astigmatism 2.09 D. Include Posterior Corneal Astigmatism checkbox is unchecked.
- Biometry:** Axial Length 24.11 mm. Method: Optical or Immersion. A-constant 119.30.
- Calculation Preferences:** SE IOL Power 20.0 D. K Index 1.3375. Refractive Cylinder Convention: Plus selected, Minus unselected.
- Final Results:** IOL Details table with columns: IOL Model, Orientation, Residual Astigmatism (Cylinder, Axis).

At the bottom of the interface, there are buttons for 'Calculate Results', 'Clear Entries', and 'Print Results'. The footer contains 'TERMS OF SERVICE' and 'LRI CALCULATOR'.

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The screenshot shows the 'Preferencias y ayuda: página de preferencias' interface. At the top, it features the TECNIS™ Family of IOLs logo and the title 'TECNIS™ Toric Calculator Europe - Middle East - Africa'. The Johnson & Johnson VISION logo is in the top right. Below the title bar, there are two tabs: 'Calculator' and 'Preferences & Support', with the latter being active. The 'Preferences' section has a header with 'Reset To Defaults' and 'Apply Preferences' buttons. It is divided into 'Display' and 'Calculation' sub-sections. Under 'Display', there are dropdowns for 'Location' (Europe - Middle East - Africa) and 'Language' (English). 'K Notation' has radio buttons for 'D' and 'mm', with 'mm' selected. 'Decimal' has radio buttons for '0.00' and '0,00', with '0.00' selected. Under 'Calculation', 'Refractive Cylinder Convention' has radio buttons for 'Plus' and 'Minus', with 'Plus' selected. 'K Index' is a dropdown menu currently set to '-Select-'. 'Method' is a dropdown menu set to 'Optical or Immersion'. There is a checkbox for 'Include Posterior Corneal Astigmatism' which is unchecked. The 'A-constant' field contains the value '119.30'. The 'Support' section below contains links for 'Training and Information Documents', including 'TECNIS™ Toric Calculator User Tutorial', 'TECNIS™ Toric Calculator Frequently Asked Questions (FAQ)', and 'TECNIS™ Toric Calculator Terms of Service'. It also includes contact information for Johnson & Johnson Surgical Vision, Inc. and AMO Ireland, along with various certification logos (MDI, CE, UDI) and a phone number.

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The screenshot shows the 'Preferencias y ayuda' page for the TECNIS Toric Calculator. The page is divided into two main sections: 'Preferences' and 'Support'. The 'Preferences' section includes options for 'Display' (Location, Language, K Notation, Decimal), 'Calculation' (Refractive Cylinder Convention, K Index), and 'Method' (Optical or Immersion, A-constant). The 'Support' section includes links for 'Training and Information Documents' and contact information for Johnson & Johnson Surgical Vision, Inc. and AMO Ireland. The page also features a footer with 'TERMS OF SERVICE' and 'LRI CALCULATOR'.

TECNIS™ Family of IOLs

TECNIS™ Toric Calculator
Europe - Middle East - Africa

Johnson & Johnson VISION

Calculator | Preferences & Support

Reset To Defaults | Apply Preferences

Preferences

Display

Location: Europe - Middle East - Africa

Language: English

K Notation: D mm

Decimal: 0.00 0,00

Calculation

Refractive Cylinder Convention: Plus Minus

K Index: -Select-

Method: Optical or Immersion

A-constant: 119.30

Include Posterior Corneal Astigmatism

Support

Training and Information Documents

[TECNIS™ Toric Calculator User Tutorial](#)

[TECNIS™ Toric Calculator Frequently Asked Questions \(FAQ\)](#)

[TECNIS™ Toric Calculator Terms of Service](#)

Johnson & Johnson Surgical Vision, Inc.
31 Technology Drive, Suite 200
Irvine, CA 92618 USA
2024-11-14

AMO Ireland
Block B
Liffey Valley Office Campus
Quarryvale, Co. Dublin, Ireland

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
Johnson & Johnson Surgical Vision, Inc.
31 Technology Drive, Suite 200,
Irvine, CA 92618 USA

AMO Ireland
Block B
Liffey Valley Office Campus
Quarryvale, Co. Dublin, Ireland

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AMO Australia Pty Ltd.
1-5 Khartoum Road,
North Ryde, NSW 2113
Australia

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© 2024 Johnson & Johnson Surgical Vision, Inc.

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AMO Switzerland GmbH
TECNIS es una marca comercial perteneciente u otorgada bajo licencia a Johnson & Johnson Surgical Vision, Inc.