# Calculateur torique TECNIS<sup>TM</sup>

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Johnson & Johnson VISION

Registration

Europe - Middle East - Africa English I have read and agree to the terms below.

Accept

TECNIS™ intraocular lenses are currently available for sale in select markets. This application is therefore intended only for residents of those markets in which the respective TECNIS™ intraocular lens(es) are available, and this application may only be used for those TECNIS™ intraocular lens(es) that are approved for use in your market. Using this mobile application and clicking the "Accept" button constitutes your acknowledgement of these restrictions, and your acknowledgement that you selected your proper location of residence.

TECNIS™ TORIC CALCULATOR LICENSE AGREEMENT, TERMS OF USE AND PRIVACY POLICY

PLEASE READ THIS USE AND LICENSE AGREEMENT ("AGREEMENT") IN FULL

TERMS OF SERVICE LRI CALCULATOR

The product-specific application is intended for the residents of a particular country or countries. As a result, the application may contain information on pharmaceuticals, medical devices and other products or uses of those products that are not approved in other countries or regions. If you are a resident of a country other than those to which the application lists, please return to surgical.jnjvision.com or contact your local Johnson & Johnson Vision affiliate to obtain the appropriate product information for your country of residence.

This application is intended for use in connection with Johnson & Johnson Vision IOLs only and should not be used in connection with any other IOLs.

©2024 Johnson & Johnson Surgical Vision, Inc.

Unless otherwise specified, all product names appearing in this program are trademarks owned by or licensed to Johnson & Johnson Vision, its subsidiaries or affiliates. No use of any Johnson & Johnson Vision trademark, trade name, or trade dress in this program may be made without the prior written authorization of Johnson & Johnson Vision, except to identify the product or service of the company.

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The screenshot displays the TECNIS Toric Calculator interface. The top left corner features the TECNIS Family of IOLs logo. The main area is divided into several sections: Surgeon and Patient Information, Keratometry, and Biometry. A central panel lists available lens options, with 'TECNIS Eyhance™ DIU' selected. To the right, a circular diagram illustrates the lens axis orientation with degree markings (0°, 45°, 90°, 135°, 180°, 225°, 270°, 315°). The bottom of the interface includes buttons for 'Calculate Results', 'Clear Entries', and 'Print Results', along with a footer for 'TERMS OF SERVICE' and 'LRI CALCULATOR'.

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The screenshot displays the TECNIS™ Toric Calculator interface. At the top, it shows the TECNIS Eyhance™ logo, the user 'DIU', the region 'Europe - Middle East - Africa', and the Johnson & Johnson VISION logo. The main interface is divided into sections: 'Surgeon and Patient Information', 'Keratometry', 'Biometry', 'Calculation Preferences', and 'Final Results'. The 'Surgeon and Patient Information' section includes fields for Surgeon Name, Date (11/14/2024), Patient Information, Patient Age, Eye Selection (OD/OS), K Notation (D/mm), and a Notes field. The 'Keratometry' section includes Surgically Induced Astigmatism (SIA), Flat K1, Steep K2, Preop Corneal Astigmatism, and an option to 'Include Posterior Corneal Astigmatism'. The 'Biometry' section includes Axial Length, Method (Optical or Immersion), and A-constant (119.30). The 'Calculation Preferences' section includes SE IOL Power, K Index, and Refractive Cylinder Convention (Plus/Minus). The 'Final Results' section includes a table for IOL Details and Residual Astigmatism. A diagram of the eye with meridians (0°, 45°, 90°, 135°, 180°, 225°, 270°, 315°) is shown on the right, along with a 'Temporally' and 'Nasally' label. The bottom of the interface has buttons for 'Calculate Results', 'Clear Entries', and 'Print Results', and a footer with 'TERMS OF SERVICE' and 'LRI CALCULATOR'.

<sup>1</sup> Koch DD et al. Contribution of posterior corneal astigmatism to total corneal astigmatism. J Cataract Refract Surg. 2012 Dec;38(12):2080-7

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The screenshot displays the TECNIS™ Toric Calculator interface. The top navigation bar includes the TECNIS Eyhance logo, the text 'DIU', and the Johnson & Johnson VISION logo. The main interface is divided into two sections: 'Calculator' and 'Preferences & Support'. The 'Calculator' section contains several input fields: 'Surgeon and Patient Information' (Surgeon Name, Date, Patient Information, Patient Age, Eye Selection), 'Keratometry' (Surgically Induced Astigmatism, Flat K1, Steep K2, Preop Corneal Astigmatism), 'Biometry' (Axial Length, Method, A-constant), and 'Calculation Preferences' (SE IOL Power, K Index, Refractive Cylinder Convention). The 'Final Results' section displays a table of IOL Model, Orientation, Cylinder, and Axis. A diagram of an eye is shown on the right, with a blue circle indicating the incision location and an orange/red arc indicating the astigmatism orientation. The diagram is labeled 'OD' and 'Nasal'.

IOL Details		Residual Astigmatism	
IOL Model	Orientation	Cylinder	Axis
<input type="radio"/> DIU300	101 °	+0.33 D	101 °
<input type="radio"/> DIU375	101 °	+0.16 D	11 °
<input type="radio"/> DIU450	101 °	+0.65 D	11 °

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DIU


[Return to Calculator](#)

Surgeon Name	Patient Information	Date	Patient Age
Surgeon	XYZ-123-ABC	2/6/2024	65

**OD (Right)**

T  
e  
m  
p  
o  
r  
a  
l

N  
a  
s  
a  
l

— (Incision Location) 180°  
— Calculated Orientation

IOL Details		Residual Refraction	
IOL Model	Orientation	Cylinder	Axis
DIU300	101 °	+0.33 D	101 °

Calculation Details		
SE IOL Power	K Index	Refractive Cylinder Convention
+20.0	1.3375	Plus

Keratometry and Biometry						
Flat K1 @ Meridian	Steep K2	SIA @ Meridian	AxL Length	Biometry Method	A-constant	PCA
41.25@12°	43.34@102°	0.30@180°	24.11	Optical or Immersion	119.30	No

IOL Details		Residual Astigmatism	
IOL Model	Orientation	Cylinder	Axis
DIU300	101 °	+0.33 D	101 °
DIU375	101 °	+0.16 D	11 °
DIU450	101 °	+0.65 D	11 °

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The screenshot shows the LRI Calculator interface with a 'Clear Entries' dialog box open. The dialog box contains the following text: 'Please confirm that you wish to clear all entries that you have made on this screen. (e.g., OD (Right) & OS (Left))'. There are 'OK' and 'Cancel' buttons at the bottom right of the dialog box.

The background interface includes the following sections:

- Surgeon and Patient:** Surgeon Name (Surgeon), Date (11/14/2024), Patient Information (XYZ-123-ABC), Patient Age (65).
- Eye Selection:** OD (Right) selected, OS (Left), K Notation (D selected, mm).
- Keratometry:** Surgically Induced Astigmatism (SIA) 0.30 D @ Meridian (Incision Location) 180°, Flat K1 41.25 D Flat K1 @ Meridian 12°, Steep K2 43.34 D Steep K2 @ Meridian 102°, Preop Corneal Astigmatism 2.09 D. Includes checkbox for 'Include Posterior Corneal Astigmatism'.
- Biometry:** Axial Length 24.11 mm, Method Optical or Immersion, A-constant 119.30.
- Calculation Preferences:** SE IOL Power 20.0 D, K Index 1.3375, Refractive Cylinder Convention Plus selected, Minus.
- Final Results:** Table with columns: IOL Model, Orientation, Residual Astigmatism (Cylinder, Axis).

At the bottom of the interface, there are buttons for 'Calculate Results', 'Clear Entries', and 'Print Results'. The footer contains 'TERMS OF SERVICE' and 'LRI CALCULATOR'.

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**TECNIS™** Toric Calculator  
Family of IOLs  
Europe - Middle East - Africa  
Johnson & Johnson VISION

Calculator | Preferences & Support

**Preferences** | Reset To Defaults | Apply Preferences

**Display**

Location: Europe - Middle East - Africa | Language: English | K Notation:  D  mm | Decimal:  0.00  0,00

**Calculation**

Refractive Cylinder Convention:  Plus  Minus | K Index: -Select- | Method: Optical or Immersion | A-constant: 119.30

Include Posterior Corneal Astigmatism

**Support**

**Training and Information Documents**

- [TECNIS™ Toric Calculator User Tutorial](#)
- [TECNIS™ Toric Calculator Frequently Asked Questions \(FAQ\)](#)
- [TECNIS™ Toric Calculator Terms of Service](#)

Johnson & Johnson Surgical Vision, Inc.  
31 Technology Drive, Suite 200  
Irvine, CA 92618 USA  
2024-11-14

AMO Ireland  
Block B  
Liffey Valley Office Campus  
Quarryvale, Co. Dublin, Ireland

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The screenshot shows the 'TECNIS™ Toric Calculator' interface. At the top, it displays 'TECNIS™ Family of IOLs' and 'Johnson & Johnson VISION'. The page is split into two tabs: 'Calculator' and 'Preferences & Support'. The 'Preferences & Support' tab is active, showing two sub-sections: 'Preferences' and 'Support'.

**Preferences Section:**

- Display:** Location (Europe - Middle East - Africa), Language (English), K Notation (D, mm), Decimal (0.00, 0,00).
- Calculation:** Refractive Cylinder Convention (Plus, Minus), Method (Optical or Immersion), K Index (-Select-), A-constant (119.30).
- Include Posterior Corneal Astigmatism


**Support Section:**

- Training and Information Documents:**
  - [TECNIS™ Toric Calculator User Tutorial](#)
  - [TECNIS™ Toric Calculator Frequently Asked Questions \(FAQ\)](#)
  - [TECNIS™ Toric Calculator Terms of Service](#)
- Contact Information:**
  - Johnson & Johnson Surgical Vision, Inc. 31 Technology Drive, Suite 200 Irvine, CA 92618 USA 2024-11-14
  - AMO Ireland Block B Liffey Valley Office Campus Quarryvale, Co. Dublin, Ireland
- Compliance:** IMD, CE 0344, UDI (01)05050474708822(8012) - 4.9.0.0
- Footer:** www.e-ifu.com EU: +353 1 7782500

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Johnson & Johnson Surgical Vision, Inc.  
31 Technology Drive, Suite 200,  
Irvine, CA 92618 USA

AMO Ireland  
Block B  
Liffey Valley Office Campus  
Quarryvale, Co. Dublin, Ireland

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AMO Australia Pty Ltd.  
1-5 Khartoum Road,  
North Ryde, NSW 2113  
Australia

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© 2024 Johnson & Johnson Surgical Vision, Inc.

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AMO Switzerland GmbH
TECNIS est une marque de commerce appartenant ou concédée sous licence à Johnson & Johnson Surgical Vision, Inc.