TECNIS™ Torik Hesaplayıcı

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Johnson & Johnson VISION

Registration

Europe - Middle East - Africa English I have read and agree to the terms below.

Accept

TECNIS™ intraocular lenses are currently available for sale in select markets. This application is therefore intended only for residents of those markets in which the respective TECNIS™ intraocular lens(es) are available, and this application may only be used for those TECNIS™ intraocular lens(es) that are approved for use in your market. Using this mobile application and clicking the "Accept" button constitutes your acknowledgement of these restrictions, and your acknowledgement that you selected your proper location of residence.

TECNIS™ TORIC CALCULATOR LICENSE AGREEMENT, TERMS OF USE AND PRIVACY POLICY

PLEASE READ THIS USE AND LICENSE AGREEMENT ("AGREEMENT") IN FULL

TERMS OF SERVICE LRI CALCULATOR

The product-specific application is intended for the residents of a particular country or countries. As a result, the application may contain information on pharmaceuticals, medical devices and other products or uses of those products that are not approved in other countries or regions. If you are a resident of a country other than those to which the application lists, please return to surgical.jnjvision.com or contact your local Johnson & Johnson Vision affiliate to obtain the appropriate product information for your country of residence.

This application is intended for use in connection with Johnson & Johnson Vision IOLs only and should not be used in connection with any other IOLs.

©2024 Johnson & Johnson Surgical Vision, Inc.

Unless otherwise specified, all product names appearing in this program are trademarks owned by or licensed to Johnson & Johnson Vision, its subsidiaries or affiliates. No use of any Johnson & Johnson Vision trademark, trade name, or trade dress in this program may be made without the prior written authorization of Johnson & Johnson Vision, except to identify the product or service of the company.

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The screenshot displays the TECNIS™ Toric Calculator interface. The top left corner features the TECNIS™ Family of IOLs logo. The main area is divided into several sections: Surgeon and Patient Information, Keratometry, Biometry, and a central lens selection list. The lens selection list is highlighted with a blue border and includes options like TECNIS Eyhance™, TECNIS Toric II 1-Piece IOL, TECNIS Toric 1-Piece IOL, TECNIS PureSee™, TECNIS Synergy™, TECNIS Symphony™, and TECNIS Multifocal. The right side of the interface shows a diagram of an eye with a lens and a circular diagram with angles (90°, 135°, 45°, 180°, 225°, 270°, 315°, 0°). The bottom right corner contains a table for IOL Details and Residual Astigmatism.

IOL Details		Residual Astigmatism	
IOL Model	Orientation	Cylinder	Axis

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TECNIS™ Toric Calculator
Europe - Middle East - Africa
CHANGE LENS

TECNIS Eyhance™ with TECNIS SIMPLICITY™ DIU Johnson & Johnson VISION

Calculator Preferences & Support

Surgeon and Patient Information (Notes...)

Surgeon Name: [] Date: 11/14/2024
Patient Information: [] Patient Age: []
Eye Selection: OD (Right) OS (Left) K Notation: D mm

Keratometry

Surgically Induced Astigmatism (SIA): [] D @ Meridian (Incision Location): []
Flat K1: [] mm Flat K1 @ Meridian: []
Steep K2: [] mm Steep K2 @ Meridian: []
Preop Corneal Astigmatism: [] D Include Posterior Corneal Astigmatism

Biometry

Axial Length: [] mm
Method: Optical or Immersion
A-constant: 119.30

Calculation Preferences

SE IOL Power: -Select- D
K Index: -Select-
Refractive Cylinder Convention: Plus Minus

Calculate Results Clear Entries Print Results

TERMS OF SERVICE LRI CALCULATOR

OD

Temporal Nasal

90° 45° 0° 315° 270° 225° 180° 135°

Final Results

IOL Details		Residual Astigmatism	
IOL Model	Orientation	Cylinder	Axis

¹ Koch DD ve ark. Contribution of posterior corneal astigmatism to total corneal astigmatism. J Cataract Refract Surg. 2012 Aralık;38(12):2080-7

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CHANGE LENS ▾

Calculator Preferences & Support

Surgeon and Patient Information Notes...

Surgeon Name: Surgeon Date: 11/14/2024

Patient Information: XYZ-123-ABC Patient Age: 65

Eye Selection: OD (Right) OS (Left) K Notation: D mm

Keratometry

Surgically Induced Astigmatism (SIA): 0.30 D @ Meridian (Incision Location): 180°

Flat K1: 41.25 D Flat K1 @ Meridian: 12°

Steep K2: 43.34 D Steep K2 @ Meridian: 102°

Preop Corneal Astigmatism: 2.09 D Include Posterior Corneal Astigmatism

Biometry

Axial Length: 24.11 mm

Method: Optical or Immersion

A-constant: 119.30

Calculation Preferences

SE IOL Power: 20.0 D

K Index: 1.3375

Refractive Cylinder Convention: Plus Minus

Calculate Results Clear Entries Print Results

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TECNIS Eyhance™ DIU Return to Calculator Johnson & Johnson VISION

Surgeon Name	Patient Information	Date	Patient Age
Surgeon	XYZ-123-ABC	2/6/2024	65

OD (Right)

— (Incision Location) 180°
— Calculated Orientation

IOL Details		Residual Refraction	
IOL Model	Orientation	Cylinder	Axis
DIU300	101 °	+0.33 D	101 °

Calculation Details		
SE IOL Power	K Index	Refractive Cylinder Convention
+20.0	1.3375	Plus

Keratometry and Biometry						
Flat K1 @ Meridian	Steep K2	SIA @ Meridian	AxL Length	Biometry Method	A-constant	PCA
41.25@12°	43.34@102°	0.30@180°	24.11	Optical or Immersion	119.30	No

IOL Details		Residual Astigmatism	
IOL Model	Orientation	Cylinder	Axis
DIU300	101 °	+0.33 D	101 °
DIU375	101 °	+0.16 D	11 °
DIU450	101 °	+0.65 D	11 °

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The screenshot displays the LRI Calculator software interface. A 'Clear Entries' dialog box is open, asking for confirmation to clear all entries. The dialog box text reads: 'Please confirm that you wish to clear all entries that you have made on this screen. (e.g., OD (Right) & OS (Left))'. There are 'OK' and 'Cancel' buttons. The background interface includes the following sections:

- Surgeon and Patient:** Surgeon Name, Surgeon, Date (11/14/2024), Patient Information (XYZ-123-ABC), Patient Age (65), Eye Selection (OD (Right) selected, OS (Left)), K Notation (D selected, mm).
- Keratometry:** Surgically Induced Astigmatism (SIA) (0.30 D, @ Meridian Location 180), Flat K1 (41.25 D, Flat K1 @ Meridian 12), Steep K2 (43.34 D, Steep K2 @ Meridian 102), Preop Corneal Astigmatism (2.09 D, Include Posterior Corneal Astigmatism).
- Biometry:** Axial Length (24.11 mm), Method (Optical or Immersion), A-constant (119.30).
- Calculation Preferences:** SEIOL Power (20.0 D), K Index (1.3375), Refractive Cylinder Convention (Plus selected, Minus).
- Final Results:** IOL Details (IOL Model, Orientation), Residual Astigmatism (Cylinder, Axis).

At the bottom of the interface, there are buttons for 'Calculate Results', 'Clear Entries', and 'Print Results'. The footer contains 'TERMS OF SERVICE' and 'LRI CALCULATOR'.

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Include Posterior Corneal Astigmatism

Support

Training and Information Documents

[TECNIS™ Toric Calculator User Tutorial](#)
[TECNIS™ Toric Calculator Frequently Asked Questions \(FAQ\)](#)
[TECNIS™ Toric Calculator Terms of Service](#)

 www.e-ifu.com EU: +353 1 7782500

Johnson & Johnson Surgical Vision, Inc.
31 Technology Drive, Suite 200
Irvine, CA 92618 USA
2024-11-14

AMO Ireland
Block B
Liffey Valley Office Campus
Quarryvale, Co. Dublin, Ireland

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
Johnson & Johnson Surgical Vision, Inc.
31 Technology Drive, Suite 200,
Irvine, CA 92618 USA

AMO Ireland
Block B
Liffey Valley Office Campus
Quarryvale, Co. Dublin, Ireland

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AMO Australia Pty Ltd.
1-5 Khartoum Road,
North Ryde, NSW 2113
Australia

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© 2024 Johnson & Johnson Surgical Vision, Inc.

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AMO Switzerland GmbH
TECNIS; Johnson & Johnson Surgical Vision, Inc. şirketine ait veya lisanslı bir ticari markadır.