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
Johnson & Johnson VISION

Registration

Europe - Middle East - Africa English I have read and agree to the terms below.

Accept

TECNIS® intraocular lenses are currently available for sale in select markets. This application is therefore intended only for residents of those markets in which the respective TECNIS® intraocular lens(es) are available and this application may only be used for those TECNIS® intraocular lens(es) that are approved for use in your market. Using this mobile application and clicking the “Accept” button constitutes your acknowledgement of these restrictions, and your acknowledgement that you selected your proper location of residence.

TECNIS® TORIC IOL CALCULATOR LICENSE AGREEMENT, TERMS OF USE AND PRIVACY POLICY

PLEASE READ THIS USE AND LICENSE AGREEMENT (“AGREEMENT”) IN FULL

PRIVACY POLICY TERMS OF SERVICE CONTACT LRI CALCULATOR

The product-specific application is intended for the residents of a particular country or countries. As a result, the application may contain information on pharmaceuticals, medical devices and other products or uses of those products that are not approved in other countries or regions. If you are a resident of a country other than those to which the application lists, please return to surgical.jnjvision.com or contact your local Johnson & Johnson Vision affiliate to obtain the appropriate product information for your country of residence.

This application is intended for use in connection with Johnson & Johnson Vision IOLs only and should not be used in connection with any other IOLs.

©2021 Johnson & Johnson Surgical Vision, Inc.

Unless otherwise specified, all product names appearing in this program are trademarks owned by or licensed to Johnson & Johnson Vision, its subsidiaries or affiliates. No use of any Johnson & Johnson Vision trademark, trade name, or trade dress in this program may be made without the prior written authorization of Johnson & Johnson Vision, except to identify the product or service of the company.

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The screenshot displays the TECNIS Family of IOLs Calculator interface. It is divided into several sections:

- Surgeon and Patient Information:** Fields for Surgeon Name, Patient Information, Eye Selection (OD/OS), and K Notation.
- Keratometry:** Fields for Surgically Induced Astigmatism (SIA), Flat K1, Steep K2, and Preop Corneal Astigmatism.
- Biometry:** Fields for Axial Length, Method (Optical or Immersion), and A-constant.
- Calculator:** A list of lens options categorized by type (MONOFOCAL, PRESBYOPIA CORRECTING) and model (Eyhance™, Synergy™, Symphony®, Multifocal).
- Visual Representation:** A diagram of an eye with a lens diagram showing meridians (0°, 45°, 90°, 135°, 180°, 225°, 270°, 315°) and the text "N a s a l".
- Results:** A table for IOL Details and Residual Astigmatism.

IOL Details		Residual Astigmatism	
IOL Model	Orientation	Cylinder	Axis

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The screenshot shows the 'TORIC CALCULATOR' interface for TECNIS Eyhance Toric II DIU. The page is titled 'United States' and has a 'CHANGE LENS' button. The main content is divided into several sections:

- Surgeon and Patient Information:** Includes fields for Surgeon Name, Date (3/24/2022), Patient Information, and Patient Age. Eye Selection is set to OD (Right), and K Notation is set to D.
- Keratometry:** Includes fields for Surgically Induced Astigmatism (SIA), Flat K1, Steep K2, and Preop Corneal Astigmatism. It also has a checkbox for 'Include Posterior Corneal Astigmatism'.
- Biometry:** Includes fields for Axial Length (mm), Method (Optical or Immersion), and A-constant (119.30).
- Calculation Preferences:** Includes dropdowns for SE IOL Power and K Index, and a radio button for Refractive Cylinder Convention (Plus/Minus).

On the right side, there is a diagram of an eye with meridians marked at 0°, 45°, 90°, 135°, 180°, 225°, 270°, and 315°. The vertical axis is labeled 'T e m p o r a l' and the horizontal axis is labeled 'N a s a l'.

At the bottom, there are buttons for 'Calculate Results', 'Clear Entries', and 'Print Results'. The footer contains links for 'PRIVACY POLICY', 'TERMS OF SERVICE', 'CONTACT', and 'LRI CALCULATOR'.

<sup>1</sup> Koch DD et al. Contribution of posterior corneal astigmatism to total corneal astigmatism. J Cataract Refract Surg. 2012 Dec;38(12):2080-7

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The screenshot displays the TECNIS Eyhance Toric II IOL Calculator interface. The top navigation bar includes the TECNIS Eyhance logo, 'Toric II DIU', 'United States', and the Johnson & Johnson VISION logo. A 'CHANGE LENS' dropdown menu is visible. The main interface is divided into several sections:

- Surgeon and Patient Information:** Includes fields for Surgeon Name (Dr. Surgeon), Date (3/24/2022), Patient Information (XYZ-123-ABC), and Patient Age (65). Eye Selection is set to OD (Right).
- Keratometry:** Shows Surgically Induced Astigmatism (SIA) as 0.35 D, Flat K1 as 41.25 D, Steep K2 as 43.34 D, and Preop Corneal Astigmatism as 2.09 D.
- Biometry:** Shows Axial Length as 24.11 mm, Method as Optical or Immersion, and A-constant as 119.30.
- Calculation Preferences:** Shows SE IOL Power as 20.0 D, K Index as 1.3375, and Refractive Cylinder Convention as Plus.
- Final Results:** A table showing IOL Model, Orientation, Residual Astigmatism (Cylinder), and Axis for three models: DIU300, DIU375, and DIU450.

On the right side, a diagram of an eye shows the orientation of the IOL (orange line) and the incision (blue line) relative to the nasal and temporal meridians. The diagram includes a 100-degree angle and various meridian markings (90°, 45°, 0°, 315°, 270°, 225°, 180°, 135°).

IOL Details		Residual Astigmatism	
IOL Model	Orientation	Cylinder	Axis
○ DIU300	100 °	+0.38 D	100 °
○ DIU375	100 °	+0.11 D	10 °
○ DIU450	100 °	+0.61 D	10 °

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Toric II  
DIU


Return to Calculator

Surgeon Name	Patient Information	Date	Patient Age
Dr. Surgeon	XYZ-123-ABC	3/24/2022	65

**OD (Right)**

T  
e  
m  
p  
o  
r  
a  
l

N  
a  
s  
a  
l

— (Incision Location) 180°  
— Calculated Orientation

IOL Details		Residual Refraction	
IOL Model	Orientation	Cylinder	Axis
DIU375	100 °	+0.11 D	10 °

Calculation Details		
SE IOL Power	K Index	Refractive Cylinder Convention
+20.0	1.3375	Plus

Keratometry and Biometry						
Flat K1 @ Meridian	Steep K2	SIA @ Meridian	AxL Length	Biometry Method	A-constant	PCA
41.25@12°	43.34@102°	0.35@180°	24.11	Optical or Immersion	119.30	No

IOL Details		Residual Astigmatism	
IOL Model	Orientation	Cylinder	Axis
DIU300	100 °	+0.38 D	100 °
DIU375	100 °	+0.11 D	10 °
DIU450	100 °	+0.61 D	10 °

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The screenshot displays the TECNIS Eyhance IOL Calculator interface. A modal dialog box titled "Clear Entries" is open, asking for confirmation to clear all data. The background interface includes the following sections:

- Surgeon and Patient Information:** Surgeon Name (Dr. Surgeon), Date (3/24/2022), Patient Information (XYZ-123-ABC), Patient Age (65).
- Eye Selection:** OD (Right) selected, OS (Left), K Notation (D selected, mm).
- Keratometry:** Surgically Induced Astigmatism (SIA) 0.35 D, Flat K1 41.25 D, Steep K2 43.34 D, Preop Corneal Astigmatism 2.09 D.
- Biometry:** Axial Length 24.11 mm, Method (Optical or Immersion), A-constant 119.30.
- Calculation Preferences:** SE IOL Power 20.0 D, K Index 1.3375, Refractive Cylinder Convention (Plus selected, Minus).
- Final Results:** Table with columns for IOL Model, Orientation, Residual Astigmatism (Cylinder, Axis).
- Diagram:** A circular diagram for Temporal and Nasal angles with markings at 90°, 135°, 45°, 180°, 225°, 270°, and 315°.

Buttons at the bottom of the interface include "Calculate Results", "Clear Entries", and "Print Results". A footer contains links for "PRIVACY POLICY", "TERMS OF SERVICE", "CONTACT", and "LRI CALCULATOR".

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Preferences Reset To Defaults Apply Preferences

**Display**

Location: United States Language: English K Notation:  D  mm  
Decimal:  0.00  0,00

**Calculation**

Refractive Cylinder Convention:  Plus  Minus K Index: -Select-  
Method: Optical or Immersion A-constant: 119.30

Include Posterior Corneal Astigmatism

**Support**

Training and Information Documents  
[TECNIS® Toric IOL Calculator User Tutorial](#)  
[TECNIS® Toric IOL Calculator FAQ](#)  
[TECNIS® Toric IOL Calculator Terms of Service](#)

(01)05050474708822(8012) - 4.5.2.0 CE

PRIVACY POLICY TERMS OF SERVICE CONTACT LRI CALCULATOR

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The screenshot displays the 'Preferences and Support' page of the TECNIS TORIC CALCULATOR. The page is divided into two main sections: 'Preferences' and 'Support'.

**Preferences Section:**

- Display:** Location is set to 'United States', Language is 'English'. K Notation is set to 'D' (radio button selected) and 'mm' (radio button unselected). Decimal is set to '0.00' (radio button selected) and '0,00' (radio button unselected).
- Calculation:** Refractive Cylinder Convention is set to 'Plus' (radio button selected) and 'Minus' (radio button unselected). K Index is set to '-Select-'. Method is set to 'Optical or Immersion'. A-constant is set to '119.30'.
- Include Posterior Corneal Astigmatism

**Support Section:**

- Training and Information Documents:**
  - [TECNIS® Toric IOL Calculator User Tutorial](#)
  - [TECNIS® Toric IOL Calculator FAQ](#)
  - [TECNIS® Toric IOL Calculator Terms of Service](#)

At the bottom right of the Support section, the version number '(01)05050474708822(8012) - 4.5.2.0' and a CE mark are displayed.

The footer of the page contains four links: [PRIVACY POLICY](#), [TERMS OF SERVICE](#), [CONTACT](#), and [LRI CALCULATOR](#).

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Johnson & Johnson Surgical Vision, Inc.  
31 Technology Drive, Suite 200,  
Irvine, CA 92618 USA

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AMO Australia Pty Ltd.  
1-5 Khartoum Road,  
North Ryde, NSW 2113  
Australia

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© Johnson & Johnson Surgical Vision, Inc., 2023.

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AMO Switzerland GmbH